



# Recombinant Human DALR anticodon-binding domain-containing protein 3 (DALRD3)

<b>Product Code</b>	CSB-BP707703HU
<b>Abbreviation</b>	DALRD3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5D0E6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MATRRLGVGE TLGALNAALG PGGPVWIKET RTRHLRSRDF LAPHRALQAR FDDGQVPEHL LHALACLQGP GVAPVLRCAP TPAGLSLQLQ RSAVFERVLS AVAAAYATPAS PASLGQRVLL HCPALRSSPC ALRLSQLRTV LVADHLARAL RAHGVCVRLV PAVRDPHMLT FLQQLRVDWP AASERASSHT LRSHALEELT SANDGRTLSP GILGRCLKE LVEEQGRTAG YDPNLDNCLV TEDLLSVLAE LQEALWHWPE DSHPGLAGAS DTGTGGCLVV HVVSCSEEFQ QQKLDLLWQK LVDKAPLRQK HLICGPVKVA GAPGTLMTAP EYYEFRHTQV CKASALKHGG DLAQDPAWTE IFGVLSVATI KFEMLSTAPQ SQLFLALADS SISTKGTKSG TFVMYNCARL ATLFESYKCS MEQGLYPTFP PVSSLDIFSLL HDEGEWLLLF NSILPFPDLL SRTAVLDCTA PGLHIAVRTE MICKFLVQLS MDFSSYNRV HILGEPRPHL FGQMFVRLQL LRAREVLHT GLAMLGLPPL SHI
<b>Source</b>	Baculovirus
<b>Target Names</b>	DALRD3
<b>Protein Names</b>	Recommended name: DALR anticodon-binding domain-containing protein 3
<b>Expression Region</b>	1-543
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.