



Recombinant Human DALR anticodon-binding domain-containing protein 3 (DALRD3)

Product Code	CSB-MP707703HU
Abbreviation	DALRD3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5D0E6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATRRLGVGE TLGALNAALG PGGPVWIKET RTRHLRSRDF LAPHRALQAR FDDGQVPEHL LHALACLQGP GVAPVLRCAP TPAGLSLQLQ RSAVFERVLS AVAAATPAS PASLGQRVLL HCPALRSSPC ALRLSQLRTV LVADHLARAL RAHGVCVRLV PAVRDPHMLT FLQQLRVDWP AASERASSHT LRSHALEELT SANDGRTLSP GILGRCLLKE LVEEQGRTAG YDPNLDNCLV TEDLLSVLAE LQEALWHWPE DSHPGLAGAS DTGTGGCLVV HVVSCSEEFQ QQKLDLLWQK LVDKAPLRQK HLICGPVKVA GAPGTLMTAP EYYEFRHTQV CKASALKHGG DLAQDPAWTE IFGVLSVATI KFEMLSTAPQ SQLFLALADS SISTKGTKSG TFVMYNCARL ATLFESYKCS MEQGLYPTFP PVSSLDIFSLL HDEGEWLLLF NSILPFPDLL SRTAVLDCTA PGLHIAVRTE MICKFLVQLS MDFSSYNRV HILGEPHRL FGQMFVRLQL LRAREVLHT GLAMLGLPPL SHI
Source	Mammalian cell
Target Names	DALRD3
Protein Names	Recommended name: DALR anticodon-binding domain-containing protein 3
Expression Region	1-543
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.