



Recombinant Chicken Transforming growth factor beta-2 (TGFB2)

Product Code	CSB-BP023448CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30371
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	LSTCSTLDMD QFMRKRIEAI RGQILSKLKL TSPPEDEYPEP EEVPPEVISI YNSTRDLLQE KANHRAATCE RERSDEEYYA KEVYKIDMQP FYPENAIPPS YYSLYFRIVR FDVSAMEKNA SNLVKAEFRV FRLQNSKARV SEQRIELYQV LKSSELSSPG QRYIDSKVVK TRAEGEWLSF DVTEAVHEWL HHRDRNLGFK ISLHCPCTF VPSNNYIIPN KSEEPPEARFA GIDDYTYSSG DVKALKSNRK KYSKGKTPHLL LMLPSYRLE SQQPSRRKKR
Source	Baculovirus
Target Names	TGFB2
Protein Names	Recommended name: Transforming growth factor beta-2 Short name= TGF-beta-2
Expression Region	21-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.