



Recombinant Staphylococcus aureus 1,4-Dihydroxy-2-naphthoyl-CoA synthase (menB)

Product Code	CSB-BP708213FLB
Abbreviation	menB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5HH38
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain COL)
Purity	>85% (SDS-PAGE)
Sequence	MTNRQWETLR EYDEIKYEFY EGIKVTINR PEVRNAFTPK TVAEMIDAFS RARDDQNVSV IVLTGEGDLA FCSGGDQKKR GHGGYVGEDQ IPRLNVLDLQ RLIRIIPKPV IAMVKGAVG GGNVLNVCD LTIAADNAIF GQTGPKVGSF DAGYGSGYLA RIVGHKKARE IWYLCRQYNA QEALDMGLVN TVVPLEKVED ETVQWCKEIM KHSPTALRFL KAAMNADTDG LAGLQQMAGD ATLLYYTTDE AKEGRDAFKE KRDPDFDQFP KFP
Source	Baculovirus
Target Names	menB
Protein Names	Recommended name: 1,4-Dihydroxy-2-naphthoyl-CoA synthase Short name= DHNA-CoA synthase EC= 4.1.3.36
Expression Region	1-273
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.