



Recombinant Human NHS-like protein 2 (NHSL2)

Product Code	CSB-MP708385HU
Abbreviation	NHSL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5HYW2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MESMGMVYSV PSSCNGPTES TFSTSWKGDA FTYMTPSATS QSNQVNENGK NPSCGNSWVS LNKVPPLVPK EAATLLVARD NPAGCSGSAG YPERLIQQRH MPERPSKIGL LTSGTSRLET GPGGASRFRE RSLSVPTDSG TTDVDYDEEQ KANEACALPF ASTSSEGSNS ADNIALSAQ QEAQHRRQRS KSISLRKAKK KPSPPTRSVS LVKDEPGLLP EGGSALPKDQ RPKSLCLSLE HQGHHSSHPD AQGHPAIPNH KDPESTQFSH HWYLTDWKSG DTYQSLSSSS TATGTTVIEC TQVQGSSESL ASPSTSRATT PSQLSIEVEA REISSPGRPP GLMSPSSGYS SQSETPTPTV SMSLTGLHLP PPSSSVRVRP VUPERKSSLP PTSPMEKFPK SRLSFDLPLT SSPNLDLSGM SISIRSKTKV SRHHSETNFG VKLAQKTNPQ QPIMPMVTQS DLRSVRLRSV SKSEPEDDIE SPEYAEPPRA EEVFTLPERK TKPPVAEKPP VARRPPSLVH KPPSVPEEYA LTSPTLAMPP RSSIQHARPL PQDSYTVVRK PKPSSFPDGR SPGESTAPSS LVFTPFASSS DAFFSGTQQP PQGSVEDEGP KVRVLPERIS LQSQEEAEKK KGKIPPPVPK KPSVLYLPLT SPTAQMEAYV AEPRLPLSPI ITLEEDTKCP ATGDDLQSLG QRVSTPQAD SEREASPLG
Source	Mammalian cell
Target Names	NHSL2
Protein Names	Recommended name: NHS-like protein 2
Expression Region	1-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.