



Recombinant *Drosophila virilis* Mediator of RNA polymerase II transcription subunit 29 (ix)

Product Code	CSB-EP708422DMP
Abbreviation	ix
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ISW5
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila virilis</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MNPNMNMMQM SGPPMMQVSP MMQSSPQPMM PTGPPGPVPM QQQHQQQQQQ QQQQQQQQQQ QAEKLDNISR VKSLLGPLRE SMFLTIRSSA FTLQQNNLAD NLKRDTGGHG HVPRFDKHL E DFYACCDQME LHLKTAIQCM QQLTSSQHYL PGAVTAMRME NFMQDNPAGP IPYPTYLNTV RVHVQSAKDI HDTLISAAQN ISQAD
Source	<i>E.coli</i>
Target Names	ix
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 29 Alternative name(s): Mediator complex subunit 29 Protein intersex
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.