



Recombinant *Oryza sativa* subsp. *japonica*

Putative protein ABIL2 (Os01g0551800, LOC_Os01g37110)

Product Code	CSB-MP708475OFG
Abbreviation	Os01g0551800, LOC_Os01g37110
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JKN2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MVEAVASFRR GGVGIAGMSS STLEEVQMQE TLIFSDTIKD LKMLKSQLYS AAEYFELAYT QEDDKQEVNM NLKEYSVKAL VNTVDHLGSI SFKVSSLIDQ RFDEVDDTNL RVSCIHQRAQ VSQACMDKEG LSQQSLVITA PKYHKRYILP AGDGSMPPNAV PNFSEMRKAK NRAAQMQQVF SAAAASQAKA KEKQPSFSKL RSIARAPSQR ARSSSPAQRP PSENTIPTKR ADKRSESP RTTPLTRSGS LPQKPSLLKT SSVRVQMHTS EHKKLASVRS QADRNDDEK EQTPKKGKKF LKSLLSRRKS RKEEPLPCYF DDY
Source	Mammalian cell
Target Names	Os01g0551800
Protein Names	Recommended name: Putative protein ABIL2 Alternative name(s): Abl interactor-like protein 2
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.