



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 12 (Os01g0917400, LOC_Os01g68860)

Product Code	CSB-BP708476OFG
Abbreviation	Os01g0917400, LOC_Os01g68860
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JLB5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDDAGRASAP AVVTVTASAA APTPPLPPPP PPPPSSLPA TAAATDEPSH DPAALYGEEM WQQMTMSGSG AMQPGPYPER SGEPCDCTYYL RTGLCRFGMS CRFNHPQDRN LAIASARMKG EYPERMGQPE CQYYLKTGTC KFGPTCKFHH PREKAGIAGR VQLNTLGYPL RPSEKECAYY LKTGQCKYGN TCKFHPELF NAMASSRGSP IYPSVHSSAT AGPPYTGTMA SWAFPRGSFI PSPRWQNPSN YAPMIVPQGL VQVPSWNSYT GQMMPVSSSE SRLQSPGAQQ TYGTSQQVDA SAGNQGMLSP YRSSYPVPQ YALQRENVFP ERPDQPECQY YMKTGDCKFG AVCKFHHPRV RSMPTDCVL SPVGLPLRPG EELCKFYSRY GICKFGANCK FDHPTMAPP MGVYAYGSAST NVPMVRLLQ SPSASAYTS
Source	Baculovirus
Target Names	Os01g0917400
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 12 Short name= OsC3H12 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 2
Expression Region	1-439
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.