



# Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 12 (Os01g0917400, LOC\_Os01g68860)

<b>Product Code</b>	CSB-MP708476OFG
<b>Abbreviation</b>	Os01g0917400, LOC_Os01g68860
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5JLB5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDDAGRASAP AVVTVTASAA APTPPLPPPP PPPPPSQLPA TAAATDEPSH DPAALYGEGM WQQMTMSGSG AMQPGPYPER SGEPDCTYYL RTGLCRFGMS CRFNHPQDRN LAIASARMKG EYPERMGQPE CQYYLKTGTC KFGPTCKFHH PREKAGIAGR VQLNTLGYPL RPSEKECAYY LKTGQCKYGN TCKFHHPPELF NAMASSRGSP IYPSVHSSAT AGPPYTGTMA SWAFPRGSFI PSPRWQNPSN YAPMIVPQGL VQVPSWNSYT GQMMPVSSSE SRLQSPGAQQ TYGTSQQVDA SAGNQGMLSP YRSSSYPVPQ YALQRENVFP ERPDQPECQY YMKTGDCKFG AVCKFHHPRV RSMPTDCVL SPVGLPLRPG EELCKFYSRY GICKFGANCK FDHPTMAPPM GVIYAYGSAST NVPMVRRLLQ SPSASAYTS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Os01g0917400
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 12 Short name= OsC3H12 Alternative name(s): Zinc finger CCCH domain-containing protein ZFN-like 2
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.