



Recombinant Human Uncharacterized protein C3orf38 (C3orf38)

Product Code	CSB-MP708479HU
Abbreviation	C3orf38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JPI3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEMSGLSFSE MEGCRNLLGL LDNDEIMALC DVTNRLVQP QDRQDAVHAI LAYSQSAEEL LRRRKVHREV IFKYLATQGI VIPPATEKHN LIQHAKDYWQ KQPQLKLNKET PEPVTKTEDI HLFQQQVKED KKAEKVDFRR LGEEFCHWFF GLLNSQNPFL GPPQDEWGPQ HFWHDVKLRF YYNTSEQNVM DYHGAEIVSL RLLSLVKEEF LFLSPNLDSH GLKCASSPHG LVMVG VAGTV HRGNTCLGIF EQIFGLIRCP FVENTWKIKF INLKIMGESS LAPGTLPKPS VKFEQSDLEA FYNVITVCGT NEVRHNVKQA SDSGTGDQV
Source	Mammalian cell
Target Names	C3orf38
Protein Names	Recommended name: Uncharacterized protein C3orf38
Expression Region	1-329
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.