



Recombinant Human Spindlin-3 (SPIN3)

Product Code	CSB-MP708489HU
Abbreviation	SPIN3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JUX0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MKTPFGKAAA GQRSRTGAGH GSVSVTMIKR KAAHKKHRSR PTSQPRGNIV GCRIQHGWDK GDEPLTQWKG TVLDQVPVNP SLYLIKYDGF DCVYGLELHR DERVSSLEVL PNRVASSRIS DTHLAEIMVG KAVEHIFETE EGSKNEWGRM VLAQAPVMNT WFYITYEKDP VLYMYQLLDD YKDGDRLRLQ DSNDSPLAER EPGEVIDSLV GKQVEYAKDD GSKRTGMVIH QVEAKPSVYF IKFDDDFHIY VYDLVKTS
Source	Mammalian cell
Target Names	SPIN3
Protein Names	Recommended name: Spindlin-3 Alternative name(s): Spindlin-like protein 3 Short name= SPIN-3
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.