



Recombinant Human EF-hand domain-containing protein 1 (EFHC1)

Product Code	CSB-MP708490HU
Abbreviation	EFHC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JVL4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVSNPVHGLP FLPGTSFKDS TKTAFHRSQT LSyrNGYAIV RRPTVGIGGD RLQFNQLSQA ELDELASKAP VLTYGQPKQA PPADFIPAHV AFDKKVLKFD AYFQEDVPMS TEEQYRIRQV NIYYLEDDS MSVIEPVVEN SGILQGKLIK RQRLAKNDRG DHYHWKDLNR GINITYGKT FRVVD CDQFT QVFLESQGIE LNPPEKMALD PYTELKQPL RKYVTPSDFD QLKQFLTFDK QVLRFYAIWD DTDSMYGECR TYIIHYYLMD DTVEIREVHE RNDGRDPFPL LMNRQRV PKV LVENAKNFPQ CVLEISDQEV LEWYTAKDFI VGKSLTILGR TFFIYDCDPF TRRYKKEKFG ITDLPRIDVS KREPPVKQE LPPYNGFGLV EDSAQNCFAL IPKAPKDV I KMLVNDNKVL RYLAVLESPI PEDKDRRFV SYFLATDMIS IFEPPVRNSG IIGGKYLGR KVVKPYSTVD NPVYYGPSDF FIGAVIEVFG HRFIILDTDE YVLKYMESNA AQYSPEALAS IQNHVRKREA PAPEAESKQT EKDPGVQELE ALIDIQKQL KDHSCDNIR EAFQIYDKEA SGYVDRDMFF KICESLNVPV DDSLVKELIR MCSHGEGKIN YYNFVRAFSN
Source	Mammalian cell
Target Names	EFHC1
Protein Names	Recommended name: EF-hand domain-containing protein 1 Alternative name(s): Myoclonin-1
Expression Region	1-640
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.