



Recombinant *Oryza sativa* subsp. *japonica* Probable protein ABIL3 (Os01g0236400, LOC_Os01g13530)

Product Code	CSB-YP708838OFG
Abbreviation	Os01g0236400, LOC_Os01g13530
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5NB83
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEAITMSPSS VSSHLDVDA ASTSEDMSSL QEGLLFSDSL KDLRNLRSQL YSAAEYFEVF YRNNSQKSTV MTSLKDYTE ALVSTVDHLG FVSYKVDNLV KERSDEVNET EFRVSSVEQR VRICQQTIDQ EGRSQQSLLI RAPKYHRRYI LPGTDIVESA IHPVSEPPRY SRQHMSRKMV KSQSISTPVG RQSTMRSARS PSPSARGTHH RSRSLSPSRK ARAKSPSPQI ISTQTKETRA GSPIPNSNPL ARSATVARRP PVHPKHFRQT SMQLHSDWSN HKEQEKSSSK GRGFLKSLLT RRRWRNDESL YSYLDEY
Source	Yeast
Target Names	Os01g0236400
Protein Names	Recommended name: Probable protein ABIL3 Alternative name(s): Abl interactor-like protein 3
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.