



Recombinant *Malus domestica* Maturase K (matK)

Product Code	CSB-BP708954MPS
Abbreviation	matK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4LAM6
Product Type	Recombinant Protein
Immunogen Species	<i>Malus domestica</i> (Apple) (<i>Pyrus malus</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEEFQGYLEL DRYQQHDFLY PLIFREIYIA LAHDHGLNRS ILLDNVGYDT KYSLLIIKRL ISRMYYQQNHL IISANDSNQN KFFGYNKNLY SQMMSEGFV IVEIPFSRRL VSSLEATEIV KSYNLRSIHS IFPFLEDKFP HLNIVSDVLI PYPIHLEILV QTLRYWVKDP SSLHLLRLFL HEYSNWNSLI TPKKIIFSKS NPRLFLLLYN SHVCEYESIL LFLRNQSSHL RLTSSGIFFE RIHFYEKKKD PVEEVFVNDP PAAILWFFKD PFMHYVRYQG KSILSSKDTP LLMNKWKYYL VNLWQCHSYV WSQPGRIYIN QLSKHSLDFL GYFSSMRPNL SVVRGQMLEN SFIMDNAMKK LDTLVPIPL IGSLAKVKFC NALGHPISKS TWADSSDFDI IDRFLHICRN LSHYYSGSSR KKSlyRIKIYI LRLSCVKTLA RKHKSTVTRTF LKRLGYKIIG RILYGRRTDS FFNLPKSFLY FEEVL
Source	Baculovirus
Target Names	matK
Protein Names	Recommended name: Maturase K Alternative name(s): Intron maturase
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.