



Recombinant Danio rerio Transcriptional adapter 2-beta (tada2b)

Product Code	CSB-MP709481DIL
Abbreviation	tada2b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q503N9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MADLGKKYCV NCLADVTLNR IRCAECQDIE LCPECFSAGA EIGNHRRWHG YQQVDGGRFS LWGPEAEGGW TSREEQSLD AIEQYGFNGW EDMAAHVGAS RTPQEVMDHY VSMYIHGNLG KACIPDSIPN RVTDHTCPSG GPLSPSLTTP LPPLDITVVE QQQLGYMPLR DDYEIEYDQE AEKLISGLSV NYDDEDIEIE MKRAHVDMYV RKLRRERQRRK NIARDYNLVP AFLGRDKKDK ERERAGGTVG VGGPGGAVGS GSGATVVPAG PLGSSTAATP KRKITKEEKG QRTKLRLALCQ FMPQREFEEF FDNMHKERM LRAKVRELQRY RRNGITRLDE SAEYEAARHK REKRKENKSI AGSKRGSSGG GGGTAGLGGG VGAGGGLGGG GGVSTIKEEG KDSEFSAIEN LSGFELLSDR EKVLCNSMNL SPMRYLTVKT IIIKDHLQKR QGIPSKSRLP SYLDKVLKKR ILNFLSESGW ISRNAS
Source	Mammalian cell
Target Names	tada2b
Protein Names	Recommended name: Transcriptional adapter 2-beta
Expression Region	1-486
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.