



Recombinant Danio rerio Charged multivesicular body protein 6 (chmp6)

Product Code	CSB-BP709482DIL
Abbreviation	chmp6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q503V0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	GNLFGKKKA TRVTEQDRAV LQLKQQRDKL KQYQKRITLQ MDKERQLAKQ LLKDGKKDKA LLLLKKKRYQ DQLLEKTENQ ISNLERMVQD IEFAQIEMKV IEGLKVGNDK LKKMHEVLSI EEEKIMDET HDAIEYQKQI DEMLAGSLTQ EDEEAVLAEL EAITQGEADL ELPEVPG EEL PEVPEQEPVR EKERVKKKPE REMVAV
Source	Baculovirus
Target Names	chmp6
Protein Names	Recommended name: Charged multivesicular body protein 6 Alternative name(s): Chromatin-modifying protein 6
Expression Region	2-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.