



Recombinant Replication initiation protein (repA)

Product Code	CSB-EP709549ENL
Abbreviation	repA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q52346
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	<p>MAGLKNTSYN AVHWSQLAPE EQIRFWEDYE AGRATTFLVE PERKRTKRRR GEHSTKPKCE NPSWYRPERY KALSGQLGHA YNRLVKKDPV TGEQSLRMHM SLHPFYVQKR TYAGRKYAFR PEKQRLLDTI WPVLVSFSDA VTHTVGM SVS RLAREISPKD SKGKVIPELE VTVSRLSRL AEQVRFVGLG VSEETLWDRE TRQRLPRYVW ITPAGWQMLG VDMVKLHEQQ QKRLRESEIR QQLIREGVLR EDEDISVHAA RKRWYLQRSQ DALKHRRAKA AASKRARRLK KLPADQQIHE MAEYLRKRLP PDEAYFCSDD HLKRLAIREL RQLELTLAAP PPH</p>
Source	E.coli
Target Names	repA
Protein Names	Recommended name: Replication initiation protein Alternative name(s): Replication-associated protein
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>