



Recombinant Dictyostelium discoideum Probable replication factor C subunit 3 (rfc3)

Product Code	CSB-MP709630DKK
Abbreviation	rfc3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54BN3
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MLWIDKYKPT SLDKMDYHND ISINLKNMIK SGDFPHLLVY GPSGAGKKTR ILAVLQEIYG PNALKLKIDH RTFKHPTSSK NIQITTISSH YHIEINPGEA GSYDRVVIQT IIKEIAQSPP IDSGSLGAFK IVILNEVDKL SKDAQHALRR TMEKYATFCR LILCCDSTAK VIDPIKSRCL GIRVPAPSQE EIEKVLAKVA TAEKFDLPSK LAVNVAKQSG GNLRYALMLL ESQKAKQYPF QSTELPLLDW ENYISQIVKD CFEEQSPARL SIVRGKLYEL LGHCIPPELI FKTLLEIFK KLDHNMKFEI IHWASYEHR SQIGSKPIFH LEAFIAKFMS VYKKYNQ
Source	Mammalian cell
Target Names	rfc3
Protein Names	Recommended name: Probable replication factor C subunit 3 Alternative name(s): Activator 1 subunit 3
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.