



# Recombinant Dictyostelium discoideum General transcription factor IIH subunit 4 (gtf2h4)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP709637DKK   |
| <b>Abbreviation</b>      | gtf2h4  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q54C29  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Dictyostelium discoideum (Slime mold)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MTLSKVFQYL ASLDSKDL EE LYKDPWTCQA ILRSLPPRSK QYILKMLLVD<br>TYPLSLAKDW STQASIQQHK ESLKKLFDLK IIFLDKINKP IQPQQQQSSQ<br>QSSSQQQQQQ QQQQQQTEQT IRLNPLFQDN IKRSLVQVNQ VIFSNNSSIK<br>DNHKPPSIDD LDSYSKSQWE KVLVFLSDDT VQPSKLISEL LLSSNLTKQE<br>GDGLSITSEG FKFLKDVYT QIWTLIVYL DDLEKKKKGK GSGSRNDLLSF<br>LFRLSFLNLG RGYLVSELS QQKEYLFALK QFGLIYMRTD SSILFYPTL<br>IISLTGKTL SLIQSISSER TQTQKEQGYI VLETNYRLYA YTSSSLQISL<br>LSLFVKMLYR LPNLAVGIIT RESIRTALIH GITADQIIDF VRHNSHPNAA<br>NSGQPIPDVV AEQILLWEAE RNRITYTKSV LYNSFPTNDC YIATLKFAKE<br>QDYIWSHDP LKTLVVKEEG NDIPIRNFIAK NFA |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | gtf2h4  |
| <b>Protein Names</b>     | Recommended name: General transcription factor IIH subunit 4 Alternative name(s): TFIH basal transcription factor complex subunit 4   |
| <b>Expression Region</b> | 1-483   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life   |



of lyophilized form is 12 months at -20°C/-80°C.