



# Recombinant Dictyostelium discoideum ESCRT-I complex subunit tsg101 (tsg101)

<b>Product Code</b>	CSB-MP709747DKK
<b>Abbreviation</b>	tsg101
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q54LJ3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Dictyostelium discoideum (Slime mold)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MYGHHGYPMH AHQQQMVNPT LAIVDNKMHT LSTFLNYIRA YRDPLRISKD LKETFHLFPN LSPFYENIPN RVNLCIKGT IPICFKGINY YLPIIVWVPL NYPQEFPTMV LDPTPEMRIV KNHHHVNLQG LVYHPYISSW SSNSTMETRV SQQQQPQPPQ NNISPPPYGS SPTNNNVAAQ QQQPPPPYGS SPSTSNSSSY TQPPPSYDSI KNKTNNTSSL PPPKQPQTSP PPPPTQQQQQ QQNNNIIPPQ QQPSPPPAYT EPTAAKQQQQ QTSPSVKPLP PLPPQPVIDK KSEMVDQCTI KLQELLSKFY DTTSLEIKDF EAHNKSLEEL SKKKLLEKEE LSNQLATYNS QIDQLNENIE QLEKWINEND KSDSNIDIDQ ILGPKDSLK QLLKLVSDDS TIEDLLYYLD KALHSNRISL EEYLKNVRSL SRDQFIIRAT VKKVQFIIRQ NQQQLQQQQQ QQNSPQRQQY NQQQYFSK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	tsg101
<b>Protein Names</b>	Recommended name: ESCRT-I complex subunit tsg101
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.