



Recombinant Dictyostelium discoideum Putative ras-related protein Rab-5B (rab5B)

Product Code	CSB-MP709817DKK
Abbreviation	rab5B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54RX9
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MSNANNNGQA TAKGGCCGGG GSQDNKANQN SGIVVSNTST DNRNKFDKGF KLVLLGEMNT GKTCISSRLV RNEFGPTDST IGA AFLVKSM VVDNINIKLE IWDTAGQERY RSLTPMYRG AAAAVVYDI TTKNTFETLK RWVSELQKHA SPNLILALAG NKVDLPNREV QVDEVNRYIS ELGNAIFFET SAASGQNINE LFTEICRRLI ENYK
Source	Mammalian cell
Target Names	rab5B
Protein Names	Recommended name: Putative ras-related protein Rab-5B
Expression Region	1-214
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.