



Recombinant *Synechococcus elongatus* Protein SrpA (srpA)

Product Code	CSB-EP709897FPY
Abbreviation	srpA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55025
Product Type	Recombinant Protein
Immunogen Species	<i>Synechococcus elongatus</i> (strain PCC 7942) (<i>Anacystis nidulans</i> R2)
Purity	>85% (SDS-PAGE)
Sequence	ATGVTPDQV LSAIEGTFGV NVGQRRNHK GTCAVGNFVA TTEAKTYSRS PLFSGQSIPV VARFSLAGGN PKAPDTAKNP RGLGLQFQLP NNRFLNMALL NTPVFGVASP EGFYENILAI RPDPTTGKPD PEKVKAFAKREK YPENKAQAAF LASNNPPTS Y ANTSYFGLHA FKFINQTNQT RLVRWQFVPQ DGEKRLTDAE LQAAPANFLE QKLIERTQDS PVKWDFWITL GQPGDAEDNP TIAWPSDRQQ VKVGTLLTA ASPQGAACE GINYDPLVLS DGIEPTNDPV LQFRSGVYAL SYSKRTRGL
Source	E.coli
Target Names	srpA
Protein Names	Recommended name: Protein SrpA EC= 1.-.-.
Expression Region	32-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.