



Recombinant Danio rerio Coenzyme Q-binding protein COQ10 homolog B, mitochondrial (coq10b)

Product Code	CSB-EP710053DIL
Abbreviation	coq10b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q567E6
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	LSSCG ILMTRAPPLF RAVSSRAVVA PARSFFNFAE SSNKRKEYSE RRIIGYSMQE MYEVVAKVED YLLFVPWCKK SDVIFRRSGF CKAKLTVGFP PVVENYTSLV STVRPHLVKA SCSDGKLFNH LETVWRFSPG LPGYPRTCTL DFAISFEFRS LLHSQLATVF FDEVVKQMVS AFERRASKLY GPETQIPREL MFHEIHHT
Source	E.coli
Target Names	coq10b
Protein Names	Recommended name: Coenzyme Q-binding protein COQ10 homolog B, mitochondrial
Expression Region	36-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.