



Recombinant Rat IST1 homolog (Ist1)

Product Code	CSB-YP710060RA
Abbreviation	Ist1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q568Z6
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MLGSGFKAER LRVNLRVIN RLKLEK KKT ELAQKARKEI ADYLAAGKDE RARIRVEHII REDYLVEAME ILELYCDLLL ARFGLIQSMK ELDSGLAESV STLIWAAPRL QSEVAELKIV ADQLCAKYSK EYGKLCRTNQ IGTVNDRMLH KLSVEAPPKI LVERYLIEIA KNYNVPYEPD SVVMAEAPVG VETDLIDVGF TDDVKKGGPG RGGGGGFTAP VGAPDGTMPM PMPMPMPMPS PSPNAPFAYP LPKGPSDFSG LPVGTYQAFP NIHPPQIPAT PPSYESVDDI NADKNVSSAQ IVGPKPEAPA KPPSRPVDNY NTFVLPELPS VPDTLPTASA GASTSASEDI DFDDLSRRFE ELK KKT
Source	Yeast
Target Names	Ist1
Protein Names	Recommended name: IST1 homolog
Expression Region	1-366
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>