



Recombinant Trypanosoma brucei brucei Flap endonuclease 1 (FEN1)

Product Code	CSB-MP710335TIY
Abbreviation	FEN1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q57WW6
Product Type	Recombinant Protein
Immunogen Species	Trypanosoma brucei brucei (strain 927/4 GUTat10.1)
Purity	>85% (SDS-PAGE)
Sequence	MGVLGLSKLL YDRTPGAIKE QELKVYFGRR IAIDASMAVY QFVIAMKGFQ EGQSVELTNE AGDVTSHLSG IFFRTLRLMID EGLRPIYVFD GKPPTLKASE LESRRQRAED AKHEFEKAKE EGDDEAMEKM SKRMVVRVGRD QMEEVKTLR LMGIPVVQAP SEAEAQCAEL VKKNKAWAVG TEDMDALAFG SRVMLRHLTY GEAKKRPIAE YHLDEILEAS GFMSQQFIDL CILLGCDYVP RISGIGPHKA WEGIKKYGSL EAFIESLDGT RYVVPEEFNY KDARNFFLEP EVTGPGEIDI QFREPDEEGL IKFLVDEKLF SKERVLKGIQ RLRDALTKKT QGRLDQFFTI TKPQKQVNSE ASTAGTKRNR GAVALPGVLQ RKSSSGHKKA VKK
Source	Mammalian cell
Target Names	FEN1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.