



Recombinant Danio rerio Protein BCCIP homolog (bccip)

Product Code	CSB-BP710666DIL
Abbreviation	bccip
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BJJ7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MASSAKRRRAV ETGQNPQDSD ESSDEGLEDS GDDSENSEEVNNEEVIVDF EAHTISDNDF HGIKTLQQL FLKCHVNTSD LTDIIQQNH IGSVIRQAEV PEDSDDEGDP DEVFGFISMV NLTERQGVLC LEKLKDMILD QCVKCSTPDG QERMENLLQG NAQSVGLLLS ERFVNVPPQI ALPLHKQLQK EMAEAQRNTK PSGKQCFC LM ISKTCKPLKK KSISAGDQAK DELLFVNDEE ECFYEQATVK FSYCVQDEAD SCATGKWSYD DEPMKPFRTV MVIPMDRMDT IMQKMTDYLT V
Source	Baculovirus
Target Names	bccip
Protein Names	Recommended name: Protein BCCIP homolog
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.