



# Recombinant Bovine Actin, gamma-enteric smooth muscle (ACTG2)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP710785BO   |
| <b>Abbreviation</b>      | ACTG2  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q5E9B5   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Bos taurus (Bovine)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | EEETALV CDNGSGLCKA GFAGDDAPRA VFPSIVGRPR HQGVMVGMGQ<br>KDSYVGDEAQ SKRGILTLKY PIEHGITNW DDMEKIWHHS FYNELRVAPE<br>EHPTLLTEAP LNPKANREKM TQIMFETFNV PAMYVAIQAV LSLYASGRRT<br>GIVLDSGDGV THNVPIYEGY ALPHAIMRLD LAGRDLTDYL MKILTERGYS<br>FVTTAEREIV RDIKEKLCYV ALDFENEMAT AASSSSLEKS YELPDGQVIT<br>IGNERFRCPE TLFQPSFIGM ESAGIHETTY NSIMKCDIDI RKDLYANNVL<br>SGGTTMYPGI ADRMQKEITA LAPSTMKIKI IAPPERKYSV WIGGSILASL<br>STFQQMWISK PEYDEAGPSI VHRKCF   |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | ACTG2  |
| <b>Protein Names</b>     | Recommended name: Actin, gamma-enteric smooth muscle<br>Alternative name(s): Gamma-2-actin Smooth muscle gamma-actin   |
| <b>Expression Region</b> | 3-376  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Target Details</b>    | Actins are highly conserved proteins that are involved in various types of cell motility, and maintenance of the cytoskeleton. In vertebrates, three main groups of actin isoforms, alpha, beta and gamma have been identified. The alpha actins are found in muscle tissues and are a major constituent of the contractile apparatus. The beta and gamma actins co-exist in most cell types as components of the cytoskeleton, and as mediators of internal cell motility. Actin, gamma 2, encoded by this gene, is a smooth muscle actin found in enteric tissues. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a   |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.