



# Recombinant Human SH3 domain-containing protein 19 (SH3D19)

<b>Product Code</b>	CSB-MP711356HU
<b>Abbreviation</b>	SH3D19
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5HYK7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNIMNTEQSQ NSIVSRIKVF EGQTNIETSG LPKKPEITPR SLPPKPTVSS GKPSVAPKPA ANRASGEWDS GTENRLKVTS KEGLTPYPPL QEAGSIPVTK PELPKKPNPG LIRSVNPEIP GRGPLAESSD SGKKVPTPAP RPLLLKKSVS SENPTYPSAP LKPVTVPPRL AGASQAKAYK SLGEGPPANP PVPVLQSKPL VDIDLISFDD DVLPTPSGNL AEESVGSEMV LDPFQLPAKT EPIKERAVQP APTRKPTVIR IPAKPGKCLH EDPQSPPLP AEKPIGNTFS TVSGKLSNVE RTRNLESNHP GQTGGFVRVP PRLPPRPVNG KTIPTQQPPT KVPPERPPPP KLSATRRSNK KLPFNRSSSD MDLQKKQSNL ATGLSKAKSQ VFKNQDPVLP PRPKPGHPLY SKYMLSVPHG IANEDIVSQN PGELSCKRGD VLVMLKQTEN NYLECQKGED TGRVHLSQMK IITPLDEHLR SRPNDPASHA QKPVDSGAPHA VVLHDFPAEQ VDDLNLTSGE IVYLLEKIDT DWYRGNCRNQ IGIFFANYVK VIIDIPEGGN GKRECVSSHC VKGSRVAVRF EYIGEQQDEL SFSEGEIILL KEYVNEEWAR GEVRGRTGIF PLNFVEPVED YPTSGANVLS TKVPLKTKKE DSGSNSQVNS LPAEWCEALH SFTAETSDDL SFKRGDRIQI LERLSDWCR GRLQDREGIF PAVFVRPCPA EAKSMLAIVP KGRKAKALYD FRGENEDEL S FKAGDIITEL ESVDLDDWMSG ELMGKSGIFP KNYIQFLQIS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	SH3D19
<b>Protein Names</b>	Recommended name: SH3 domain-containing protein 19 Alternative name(s): ADAM-binding protein Eve-1 EEN-binding protein Short name= EBP
<b>Expression Region</b>	1-790
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.