



Recombinant Pan troglodytes Beta-defensin 106A (DEFB106A)

Product Code	CSB-MP711382EQV
Abbreviation	DEFB106A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5IAB3
Product Type	Recombinant Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Purity	>85% (SDS-PAGE)
Sequence	FFDEKCNKLNK GTCKNNGGKN EELIALCQKS LKCCRTIQPC GSIID
Source	Mammalian cell
Target Names	DEFB106A
Protein Names	Recommended name: Beta-defensin 106A Alternative name(s): Beta-defensin 6 Short name= BD-6 Short name= DEFB-6 Short name= cBD-6 Defensin, beta 106 Defensin, beta 106A
Expression Region	21-65
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.