



Recombinant *Xenopus tropicalis* Serpin B5 (serpinb5)

Product Code	CSB-MP711726XBF
Abbreviation	serpinb5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5M8J5
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MDALRLANTA LAVDIFKKLC EKSATDNFVC SPLCISSSL S LIRKGSQGNT ASELEKALHF EKVKDPDFGF QLLSSDISKI SSANSLKLLK RVYVDNSIEC KKDFINSAKK PYPLELETID FKSQAEEART QINSSVKELT DGNFETVLNE GSCDENTKII MLGAASFKGK WVYTFNKSET KEMDFHINKK ETKPVQMMHL EARLSIGYIN ELKTMVLEMP FQSKHF SMLI LLPKDIEDDS TGLKKLEQDM TFEKYTHWTN PSMMAN SKVK VSLPKFKMEN SYDLKDMLKS LGINDAFNEE ASDFSEM TES KGISISQAIQ KACIEVDEDG TESADVSMER RLMNKEEFLA DHPFIYILRH NKTRTIIMLG RYCGPSEAS
Source	Mammalian cell
Target Names	serpinb5
Protein Names	Recommended name: Serpin B5
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.