



Recombinant Rat Caprin-1 (Caprin1)

Product Code	CSB-YP711733RA
Abbreviation	Caprin1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5M9G3
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	PSATSHSGS GSKSSGPPPP SGSSGSEAAA GAAAPASQHP ATGTGAVQTE AMKQILGVID KKLRLNLEKKK GKLEDDYQERM NKGRLNLDQD LDAVSKYQEV TNNLEFAKEL QRSFMALSQD IQKTIKKTAR REQLMREEAE QKRLKTVLEL QYVLDKLGDD DVRTDLKQGL SGVPILSEEE LSLLEDFYKL VDPERDMSLR LNEQYEHASI HLWDLLEGKE KPVCCTTYKA LKEIVERVFQ SNYFDSTHNS QNGLCEEEEA ASAPIVEDQV AEAPEPEPTEE YTEQSEVEST EYVNRQFMAE TQFSSGEKEQ VDEWAVETVE VVNSLQQQPQ AASPSVPEPH SLTPVAQSDP LVRRQRVQDL MAQMCGPYNF IQDSMLDFEN QTLDPAlVSA QPMNPTQNMD MPQLVCPQVH SESRLAASNQ VPVQPEATQV PLVSSTSEGY TASQPLYQPS HAAEQRPQKE PIDQIQATIS LNTEQTTASS SLPAASQPQV FQAGASKPLH SSGINVNAAP FQSMQTVFNM NAPVPPVNEP ETLKQQSQYQ ASYNQSFSSQ PHQVEQTELQ QDQLQTVVGT YHGSQDQPHQ VPGNHQQPPQ QSTGFPRSSQ PYYNSRGVSR GGSRGARGLM NGRYPANGF RGGYDGYRSS FSNTPNNGYT QSQFNAPRDY SGYQRDGYQQ NFKRGSQSG PRGAPRGRGG PPRPNRGMPQ MNTQQVN
Source	Yeast
Target Names	Caprin1
Protein Names	Recommended name: Caprin-1 Alternative name(s): Cytoplasmic activation- and proliferation-associated protein 1 GPI-anchored protein p137 Short name= GPI-p137 Short name= p137GPI RNA granule protein 105
Expression Region	2-707
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.