



Recombinant *Oryza sativa* subsp. japonica Probable indole-3-acetic acid-amido synthetase GH3.3 (GH3.3)

Product Code	CSB-BP711786OFG
Abbreviation	GH3.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5NAZ7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MLEKKATRST RVDGVSGEAV IEEFERVTRD AANVQRETLR RILAENGGVE YLRGLGLAGA TDPATFRARV PLATHADLEP YIDRIADGDA SPVLTAKPAT SISLSSGTTQ GKRKYLLFNE ELVKSTMQIY RISYAFRNRE FPVENGGKALQ FIYSSRETRT KGGLTATTAT TNVYRSEEFK ATMRDIQSQC CSPDEVIFGP DFAQSLYCHL LAGLLAAGDV QIVSATFAHS VVLAFTFER AWEDLCADIR RGEVSPSRVT SPAVRRAMAA LLAAPNPGLA DEVARKCAAL SNWYGVIPAL WPNARYVYGI MTGSMEHYVK KLRHYAGGLP LVAAEYGASE GWVGANVEPG TPPERATFTV LPDIAYFEFI PLKPVAGDGG YAEAEPVGLT EVAAGELYEV VMTTFAGNTR SSSSCMTLVA YYYYLQSKKWM NICRFCISVH ETSRNLVTCT FAATGHIMFA QK
Source	Baculovirus
Target Names	GH3.3
Protein Names	Recommended name: Probable indole-3-acetic acid-amido synthetase GH3.3 EC= 6.3.2.- Alternative name(s): Auxin-responsive GH3-like protein 3 Short name= OsGH3-3
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.