



Recombinant Anaplasma marginale Lipoyl synthase (lipA)

| | |
|--------------------------|---|
| Product Code | CSB-EP712503AAAC |
| Abbreviation | lipA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q5PAD4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Anaplasma marginale (strain St. Maries) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MENKPDWLRV RMPGGAVFDE VTELVRRYEL NTVCEEAACP NIGECWNKRH ATIMIMGSVC TRACAFCNVK VGVPNALDPD EPERVGAAIS KLGLKHVVIT SVDRDDLDPG GASHFAKCVR EIRKRSSVT VEILTPDFLG KPCAIDIIAS ARPDVYNHNL ETVPRLYSRV RPRAKYFNSL NLLKEVKDKS PGVFTKSGFM LGLGESKEEV YQVMDDLRCAGVDFIVMGQY LQPTKNNIEV HRYVPPSEFE QYKLMAYAKG FSMVAASPLA RSSYHAEDDF RKLKELRFAR AKVNE |
| Source | E.coli |
| Target Names | lipA |
| Protein Names | Recommended name: Lipoyl synthase EC= 2.8.1.8 Alternative name(s): Lip-syn Short name= LS Lipoate synthase Lipoic acid synthase Sulfur insertion protein LipA |
| Expression Region | 1-295 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |