



# Recombinant Human Breast cancer metastasis-suppressor 1-like protein (BRMS1L)

<b>Product Code</b>	CSB-BP712727HU
<b>Abbreviation</b>	BRMS1L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5PSV4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPVHSRGDKK ETNHHDEMEV DYAENEGSSS EDEDESSSV SEDGDSSEMD DEDCERRRME CLDEMSNLEK QFTDLKDQLY KERLSQVDAK LQEVIAGKAP EYLEPLATLQ ENMQIRTKVA GIYRELCLES VKNKYECEIQ ASRQHCESEK LLLYDTVQSE LEEKIRRLLEE DRHSIDITSE LWNDELQSRK KRKDPFSPDK KKPVVVSGPY IVYMLQDLDI LEDWTTIRKA MATLGPHRVK TEPPVKLEKH LHSARSEEGR LYYDGEWYIR GQTICIDKKD ECPTSAVITT INHDEVWFKR PDGSKSKLYI SQLQKGKYSI KHS
<b>Source</b>	Baculovirus
<b>Target Names</b>	BRMS1L
<b>Protein Names</b>	Recommended name: Breast cancer metastasis-suppressor 1-like protein Alternative name(s): BRMS1-homolog protein p40 BRMS1-like protein p40
<b>Expression Region</b>	1-323
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein shows sequence similarity to the human breast carcinoma metastasis suppressor (BRMS1) protein and the mammalian Sds3 (suppressor of defective silencing 3) proteins. This protein is a component of the mSin3a family of histone deacetylase complexes (HDAC).
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.