



# Recombinant Pan troglodytes Alcohol dehydrogenase 1B (ADH1B)

<b>Product Code</b>	CSB-EP712800EQV-B
<b>Abbreviation</b>	ADH1B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5R1W2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pan troglodytes (Chimpanzee)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	STAGKVIKC KAAVLWEVKK PFSIEDVEVA PPKAYEVRİK MVAVGICRTD DHHVSGNLVT PLPAILGHEA AGIVESVGEG VTTVKPGDKV IPLFTPQCGK CRVCKNPESN YCLKNDLGNP RGTLQDGTRR FTCRGKPIHH FLGTSTFSQY TVVDENAVAK IDAASPLEKV CLIGCGFSTG YGSAVNVAKV TPGSTCAVFG LGGVGLSAVM GCKAAGAARI IAVDINKDKF AKAKELGATE CINPQDYKKP IQEVLKEMTD GGVDIFSFEVI GRLDTMMASL LCCHEACGTS VIVGVPPASQ NLSINPMLLL TGRTWKGAVY GGFKSKEGIP KLVADFMakk FSLDALITHV LPFEKINEGF DLLHSGKSIR TVLTF
<b>Source</b>	E.coli
<b>Target Names</b>	ADH1B
<b>Protein Names</b>	Recommended name: Alcohol dehydrogenase 1B EC= 1.1.1.1 Alternative name(s): Alcohol dehydrogenase subunit beta
<b>Expression Region</b>	2-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.