



Recombinant Pongo abelii STE20-related kinase adapter protein alpha (STRADA)

Product Code	CSB-EP712952PYX
Abbreviation	STRADA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RBJ6
Product Type	Recombinant Protein
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Purity	>85% (SDS-PAGE)
Sequence	MSFLVSKPER IRRWSEKFI VEGLRDLELF GEQPPGDTRR KTNDA SSES I APFSKQEV MG SFLPEGGCYE LLTVIGKGF E DLMTVNLARY KPTGEYVTVR RINLEACSNE MVTFLQGELH VSKLFNHPNI VPYRATFIAD NELWVVT SFM AYGSAKDLIC THFMDGMNEL AIAYILQGV L KALDYIHHMG YVHRSVKASH ILISVDGKVY LSGLR SNLSM ISHGQRQRVV HDFPKYSIKV LPWLSPEVLQ QNLQGYDAKS DIYSVGITAC ELANGHV PFK DMPATQTLLE KLNGTVPCLL DTSTIPAEEL TMSPSRSVAN SGLSDSLTTS TPRPSNGDSP SHPYHRTFSP HFHHFVEQCL QRNP DARPSA STLLNHSFFK QIKRRASEAL PELLRPVTPI TNFEQSQSQD HSGIFGLVTN LEELEVDDWE F
Source	E.coli
Target Names	STRADA
Protein Names	Recommended name: STE20-related kinase adapter protein alpha Short name= STRAD alpha Alternative name(s): STE20-related adapter protein
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.