



# Recombinant Pongo abelii 14-3-3 protein gamma (YWHAG)

<b>Product Code</b>	CSB-EP712963PYX-B
<b>Abbreviation</b>	YWHAG
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5RC20
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVDREQLVQK ARLAEQAERY DDMAAAMKNV TELNEPLSNE ERNLLSVAYK NVVGARRSSW RVISSIEQKT SADGNEKKIE MVRAYREKIE KELEAVCQDV LSLLDNYLIK NCSETQYESK VFYLMKMGDY YRYLAEVATG EKRA TVV E S S EKAYSEAHEI SKEHMQP T H P I R L G L A L N Y S V F Y Y E I Q N A P E Q A C H L A K T A FDDAIAELDT LNEDSYKDST LIMQLLRDNL TLWTSDQQDD DGGEGNN
<b>Source</b>	E.coli
<b>Target Names</b>	YWHAG
<b>Protein Names</b>	Recommended name: 14-3-3 protein gamma Cleaved into the following chain: 1. 14-3-3 protein gamma, N-terminally processed
<b>Expression Region</b>	1-247
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.