



Recombinant Mouse Recombining binding protein suppressor of hairless (Rbpj)

Product Code	CSB-BP019486MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31266
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MPSGFQSPR TSPRARPKTR ITGALPMDYS EGLSAEERPA HAPSAGKFGE RPPPKRLTRE AMRNYLKERG DQTVLILHAK VAQKSYGNEK RFFCPPPCVY LMGSGWKKKK EQMERDGCSE QESQPCAFIG IGNSDQEMQQ LNLEGKNYCT AKTLYISDSD KRKHFMLSVK MFYGNSSDDIG VFLSKRIKVI SKPSSKKKQSL KNADLCIASG TKVALFNRLR SQTVSTRYLH VEGGNFHASS QQWGAFYIHL LDDDESEGEE FTVRDGYIHY GQTVKLVCSV TGMALPRLII RKVDKQTALL DADDPVSQLH KCAFYLKDETE RMYLCLSQER IIQFQATPCP KEQNKEMIND GASWTIISTD KAEYTFYEGM GPVLAPVTPV PVVESLQLNG GGDVAMLELT GQNFTPnlrv WFGDVEAETM YRCGESMLCV VPDISAFREG WRWVRQPvqv PVTlvrndgv IYSTSLTFTY TPEPGPRPHC SAAGAILRAN SSQVPSNESN TNSEGNyTNA STNSTSVTSS TATVVS</p>
Source	Baculovirus
Target Names	Rbpj
Protein Names	Recommended name: Recombining binding protein suppressor of hairless Alternative name(s): J kappa-recombination signal-binding protein RBP-J kappa
Expression Region	1-526
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.