



Recombinant Human Protein FAM132A (C1QTNF12)

Product Code	CSB-YP713130HU
Abbreviation	FAM132A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T7M4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GGVGARREQ RTQQPGQRAD PPNATASASS REGLPEAPKP SQASGPEFSD AHMTWLNFRV RPDDGALRKR CGSRDCKPRD LFGPPGPPGA EVTAETLLHE FQELLKEATE RRFSGLLDPL LPQGAGLRLV GEAFHCRLQG PRRVDKRTL V ELHGFQAPAA QGAFLRGSGL SLASGRFTAP VSGIFQFSAS LHVHSELQG KARLRARDVV CVLICIESLC QRHTCLEAVS GLESNSRVFT LQVQGLLQLQ AGQYASVFVD NGSGAVLTIQ AGSSFSGLLL GT
Source	Yeast
Target Names	C1QTNF12
Protein Names	Recommended name: Protein FAM132A
Expression Region	21-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.