



Recombinant Human RAB6-interacting golgin (GORAB)

Product Code	CSB-MP713131HU
Abbreviation	GORAB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5T7V8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSWAAVLAVA AARFGHFWGC RWPGPMAQGW AGFSEELRR LKQTKDPFEP QRRLPAKKSQ QQLQREKALV EQSQKLGQLD GSTSLLPEQL LSAPKQRVNV QKPPFSSPTL PSHFTLTSPV GDGQPQGIES QPKELGLENS HDGHNNVEIL PPKPDCKLEK KKVELQEKSQ WEVLQQEQRL MEEKNKRKKA LLAKAIAERS KRTQAETMKL KRIQKELQAL DDMVSADIGI LRNRIDQASL DYSYARKRFD RAEAEYIAAK LDIQRKTEIK EQLTEHLCTI IQQNELRKAK KLEELMQQLD VEADEETLEL EVEVERLLHE QEVESRRPVV RLERPFQPAE ESVTLEFAKE NRKCQEQAVS PKVDDQCGNS SSIPFLSPNC PNQEGNDISA ALAT
Source	Mammalian cell
Target Names	GORAB
Protein Names	Recommended name: RAB6-interacting golgin Alternative name(s): N-terminal kinase-like-binding protein 1 Short name= NTKL-BP1 Short name= NTKL-binding protein 1 Short name= hNTKL-BP1 SCY1-like 1-binding protein 1 Short
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.