



Recombinant Danio rerio PCI domain-containing protein 2 (pcid2)

Product Code	CSB-MP713166DIL
Abbreviation	pcid2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U3P0
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAHITLNQYL QQVLEAIDSR DGSFCAEMMS FKHPHVANPR LQLSSPEEKC QQLLEPPYDE MVA AHLRCTF AVSNHDFVEA YKCQTVVVQS FLKAFQAHKE ENWALPIMFA VTLDLRIFAN NAEQHLQQKG KGKVGDMLEK AAEQLMSCFR VCASDNRAGI DDSKKWGMLF LINQLFKIYF KINKLHLCKP LIRAIDSSNL KDEYTMAQRV TYKYYVGRKA MFDSDFKPAE ECLSFSFTHC HRSCQRNKRL ILIYLLPVKM LLGHMPHQL LKKYDLMQFA DVTRAVSEGN LLLLNAALVK HETFFIRCGI FLILEKLIKII TYRNLFKKVY HLLRTHQLPL AAFLVSLQMT KVEDVDIDEV QCILANLIYM GHIKGYISHQ HQKLVVSKQN PFPPLSSIS
Source	Mammalian cell
Target Names	pcid2
Protein Names	Recommended name: PCI domain-containing protein 2 Alternative name(s): CSN12-like protein
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.