



Recombinant Human Activating transcription factor 7-interacting protein 2 (ATF7IP2)

Product Code	CSB-MP713182HU
Abbreviation	ATF7IP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5U623
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASPDRSKRK ILKAKKT MPL SCRKQVEMLN KSRNVEALKT AIGSNVPSGN QSFSPSVITR TTEITKCSPS ENGASSLDSN KNSISEKSKV FSQNCIKPVE EIVHSETKLE QVVC SYQKPS RTTESPSRVF TEEAKDSLNT SENDSEHQTN VTRSLFEHEG ACSLKSSCCP PSVLSGVVQM PESTVTSTVG DKKTDQMVFH LETNSNSESH DKRQSDNILC SEDSGFVPVE KTPNLVNSVT SNNCADDILK TDECSRTSIS NCESADSTWQ SSLDTNNSH YQKKRMFSEN EENVKRMKTS EQINENICVS LERQTAFLEQ VRHLIQQEIY SINYELFDKK LKELNQRIGK TECRNKHEGI ADKLLAKIAK LQRRIKTVLL FQRNCLKPNM LSSNGASKVA NSEAMILDKN LESVNSPIEK SSVNYEPSNP SEKGSKKINL SSDQNKSVSE SNDDVMLIS VESPNLTPPI TSNPTDTRKI TSGNSSNSPN AEVMAVQKKL DSIIDLTKEG LSNCNTESPV SPLESHSKAA SNSKETTPLA QNAVQVPESF EHLPLPEPP APLPELVDKT RDTLPPQKPE LKVKR VFRPN GIALTWNITK INPKCAPVES YHLFLCHENS NNKLIWKKIG EIKALPLPMA CTLSQFLASN RYYFTVQSKD IFGRYGPFC D IKSIPGFSEN LT
Source	Mammalian cell
Target Names	ATF7IP2
Protein Names	Recommended name: Activating transcription factor 7-interacting protein 2 Short name= ATF7-interacting protein 2 Alternative name(s): MBD1-containing chromatin-associated factor 2
Expression Region	1-682
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.