



Recombinant Human Tetratricopeptide repeat protein 39B (TTC39B)

Product Code	CSB-BP713348HU
Abbreviation	TTC39B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VTQ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDAVLACRLR GRGNRVAALR PRPRPGGSAG PSPFALLCAG LSPEPRAGVG SEFFAWFLGG SSQRRNMALL GSRAELEADE DVFEDALETI SISSHSDMAT SSLHFASCDT QQAPRQRGAS TVSSSSSTKV DLKSGLEECA VALNFLSNK FTDALELLRP WAKESMYHAL GYSTIVVLQA VLTFEQQDIQ NGISAMKDAL QTCQKYRKKY TVVESFSSLL SRGSLEQLSE EEMHAEICYA ECLLQKAALT FVQDENMINF IKGGLKIRTS YQIYKECLSI LHEIQKNKIQ QEFFYEFEGG VKLGSGAFNL MLSLLPARI RLLLEFIGFSG NRELGLLQLR EGASGRSMRS ALCCLTILAF HTYISLILGT GEVNVAEAER LLAPFLQQFP NGSLVLFYHA RIELLKGNLE EAQEVFQKCI SVQEEWKQFH HLCYWELMWI NVFQQNWMQA YYYSDDLCKE SKWSKATYVF LKAAILSMLP EEDVVATNEN VVTLFRQVDS LKQRIAGKSI PTEKFAVRKA RRYASALPAP VKLILPALEM MYVWNGFSIV SKRKDLSEN LVTVEKAEAA LQSQNFNSFS VDDECLVKLL KGCCLKNLQR PLQAELCYNH VVESEKLLKY DHYLVPTLF ELASLYKSQG EIDKAIFLE TARNNYKDYS LESRLHFRIQ AALHLWRKPS SD
Source	Baculovirus
Target Names	TTC39B
Protein Names	Recommended name: Tetratricopeptide repeat protein 39B Short name= TPR repeat protein 39B
Expression Region	1-682
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.