



# Recombinant Human Tetratricopeptide repeat protein 39B (TTC39B)

<b>Product Code</b>	CSB-EP713348HU
<b>Abbreviation</b>	TTC39B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5VTQ0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDAVLACRLR GRGNRVAALR PRPRPGGSAG PSPFALLCAG LSPEPRAGVG SEFFAWFLGG SSQRRNMALL GSRAELEADE DVFEDALETI SISSHSDMAT SSLHFASCDT QQAPRQRGAS TVSSSSSTKV DLKSGLEECA VALNLFSLNK FTDALELLRP WAKESMYHAL GYSTIVVLQA VLTFEQQDIQ NGISAMKDAL QTCQKYRKKY TVVESFSSLL SRGSLEQLSE EEMHAEICYA ECLLQKAALT FVQDENMINF IKGGLKIRTS YQIYKECLSI LHEIQKNKIQ QEFFYEFEGG VKLGSGAFNL MLSLLPARI RLEFIGFSG NRELGLLQLR EGASGRSMRS ALCCLTILAF HTYISLILGT GEVNVAEAER LLAPFLQQFP NGSLVLFYHA RIELLKGNLE EAQEVFQKCI SVQEEWKQFH HLCYWELMWI NVFQQNWMQA YYSDLLCKE SKWSKATYVF LKAAILSMLP EEDVVATNEN VVTLFRQVDS LKQRIAGKSI PTEKFAVRKA RRYASALPAP VKLILPALEM MYVWNGFSIV SKRKDLSNL LVTVEKAEAA LQSQNFNSFS VDDECLVKLL KGCCLKNLQR PLQAELCYNH VVESEKLLKY DHYLVPTLF ELASLYKSQG EIDKAIFLE TARNNYKDYS LESRLHFRIQ AALHLWRKPS SD
<b>Source</b>	E.coli
<b>Target Names</b>	TTC39B
<b>Protein Names</b>	Recommended name: Tetratricopeptide repeat protein 39B Short name= TPR repeat protein 39B
<b>Expression Region</b>	1-682
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.