



Recombinant Human Sorting nexin-30 (SNX30)

Product Code	CSB-EP713353HU-B
Abbreviation	SNX30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5VWJ9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGGPPKALP STGPHSLRDM PHPLAGSSSE EAVGGDSTPS PDLLMARSFG DKDLILPNGG TPAGTSSPAS SSSLLNRLQL DDDIDGETRD Lfvivddpkk HVCTMETYIT YRITTKSTRV EFDLPEYSVR RRYQDFDWLR SKLEESQPTH LIPPLPEKFV VKGVVDRFSE EFVETRRKAL DKFLKRITDH PVLSFNEHFN IFLTAKDLNA YKKQGIALLT RMGESVKHVT GGYKLRTRPL EFAAIGDYLD TFALKLGTID RIAQRIKEE IEYLVELREY GPVYSTWSAL EGELAEPLG VSACIGNCST ALEELTDDMT EDPLPVLREY ILYSDSMKSV LKKRDQVQAE YEAKLEAVAL RKEDRPKVPA DVEKCQDRME CFNADLKADM ERWQNNKRQD FRQLLMGMAD KNIQYYEKCL MAWESIPLL QEKQEAQ
Source	E.coli
Target Names	SNX30
Protein Names	Recommended name: Sorting nexin-30
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.