



Recombinant *Oryza sativa* subsp. japonica E3 ubiquitin-protein ligase makorin (MKRN)

| | |
|--------------------------|---|
| Product Code | CSB-EP713694OFG |
| Abbreviation | MKRN |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q5ZA07 |
| Product Type | Recombinant Protein |
| Immunogen Species | <i>Oryza sativa</i> subsp. japonica (Rice) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSTKRVLCKF FMHGACLKGE YCEFSHDWND QPNNVCTFYQ KGSCSYGSRC RYDHVKVSRN PTVAPPPSSS TTTRASSSLQ PLSFGRPHHV GYQADSSNPR QQISMDVLAH SGSKPVWRND FQHESVLEDG IDWSISPTVQ NQTTLSPADL PICSFAAGGN CPYGEECPQM HGDLCTTCGK MCLHPYRPDE REEHTKLCEK NHRLESLKR SQEIECSVCL DRVLSKPTAA ERKFGLLSEC DHPFCISCIR NWRNNSPTSG MDVNSALRAC PICRKLSSYYV IPSVLWYFSK EEKLEIIDNY KAKLKSIDCK YDFDGTGTCP FGSSCFYKHA YRDGRLEEVI LRHLDADDGS TVIAKNIRLS DFLSRLHL |
| Source | E.coli |
| Target Names | MKRN |
| Protein Names | Recommended name: E3 ubiquitin-protein ligase makorin EC= 6.3.2.- Alternative name(s): Zinc finger CCCH domain-containing protein 41 Short name= OsC3H41 |
| Expression Region | 1-368 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.