



Recombinant Chicken THAP domain-containing protein 5 (THAP5)

Product Code	CSB-MP713697CH
Abbreviation	THAP5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZHN5
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MPRYCAASYC KNRGGQSARD QRKLSFYPPF LHDKERLEKW LRNMKRDAWT PSKHQLLCSD HFTPDSLQVW WGIRYLKHTA VPTIFSSPDD EEKGSSQNSP QEIRREDQEE TTKNVESKKV PVSLELCTPK KSSVTAENPD KKAIEVICSTV LSKSLQVQKL QLENREDFEA DSIILDNSSQ QHIHQPEPVL MAAAVQNVEA TDVHASVEVP VSCTATVLQF SDPDYLNSSL KLKNTLGSIT DYTLENPNSH VVGCSVEVQP ATENAVLVNA VTQTIEQLSG SEESVIAIIV PAESPKEPEI INRSFLPIKQ EFLDTEEAET DNSLHMNAYN GLEILQTEHS YCKQDIDREH LWQKISKLHS KITLLEMQEV KTLGRLRSLE ALIGQLKQEN LLSEEKLVKIV ENCFTTLEVT MIQ
Source	Mammalian cell
Target Names	THAP5
Protein Names	Recommended name: THAP domain-containing protein 5
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.