



Recombinant Chicken Protein FAM49A (FAM49A)

Product Code	CSB-EP713700CH-B
Abbreviation	FAM49A
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5ZI04
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MGNLLKVLTR EIENYPHFFL DFENAQPTDG EREVWNQISA VLQDSESMILA DLQAYKGAGQ EIRDAIQNPN DIQLQEKAWN SVCPLVVRLK RFYEFSLRLE KALQSLLESL TCPPTYPTQH LEREQALAKE FAELHFTLR FDELKMRNPA IQNDFSYYRR TISRNRINNM HLDIENEVNN EMANRMSLFY AEATPMLKTL SNATTHFVSE NKTLPIENTT DCLSTMASVC KVMLETPEYR SRFTSEETLM FCMRVMVGVI ILYDHVHPVG AFSKTSKIDM KGCIKVLKEQ PPDTVEGLLN ALRFTTKHLN DESTSKQIRA MLQ
Source	E.coli
Target Names	FAM49A
Protein Names	Recommended name: Protein FAM49A
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>