



Recombinant Chicken Nucleoporin-like protein 2 (NUPL2)

Product Code	CSB-BP713701CH
Abbreviation	NUPL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZI22
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MTICQFFLQG RCRFGDRCWN EHPRGGGGRP HSAGPVRGAG GGWGAASQRY ANVIQPIFK HSTWGGSGDG GGFSGASDFG SPGNKSAVFS QNRFSALSSA HPADGFSDEE QRLLDCVAKD MATWESSGQW MFSCYSPEAG KPNVSGFREF SAEVRLEY NCSANNNTEN YINSVNQLVQ ERRNRLQELK ALNASGKESL LSQKNAV TQPLPSLGFGGQ QASSFGFPSF PVSSSSGAAS FSFKANPSVP PGNAAAVGSSAAASNPTFG VTSSPSVPNP VGSGNSSAPS AASFSGKTSG TTSGCGTSGLSGFGSSAAAN SSSTAPLPVS ATPSAATGTS QSGASSASAA QTAGASGHNVTSAPSAVPNG IASDKLYTPR SELTAELEQ FEAKRFTLGK IPLKPPIDL LYL
Source	Baculovirus
Target Names	NUPL2
Protein Names	Recommended name: Nucleoporin-like protein 2
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.