



Recombinant Chicken UAP56-interacting factor (FYTTD1)

Product Code	CSB-EP713718CH
Abbreviation	FYTTD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZJ20
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MSGFGAAALL SGSSAAAGTR SGSSDSLEKI DMSLDDIIKL NKKEERKQYS PKMKRGLQQN RTQQFRTPGS KWGIQQQKGY GKNHLGHRKK IVGKKRPYGV ITGLAAKKAM GSHKGISPLN RQPLSEKNTQ RNYPIKPKP NLQRQTEMQR KQIPALRRPA PLSRRIGNKL NQQKDTRQAT FLFRRGLKVQ AQVQPTEDLD NQAAKRTRQW RTSTTSGGIL TVSIENPGAI ITPISQKLRL TRTPVPPFLM KKDQSEEKI PKGVPLQFDI NSVGKQTGMT LNERFGILKE QRTALSQNKG SRFVTVG
Source	E.coli
Target Names	FYTTD1
Protein Names	Recommended name: UAP56-interacting factor Alternative name(s): Forty-two-three domain-containing protein 1 Short name= Protein 40-2-3
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.